



# Clonazepam

## 0.1mg/mL Oral Suspension

Last updated: April 2007

Batch No: \_\_\_\_\_

Disclaimers: CAUTION: The information on The Hospital for Sick Children World Wide Web Site has been developed specifically for The Hospital for Sick Children. It is provided for informational purposes for qualified health care professionals. This material is not intended as a substitute for consulting qualified health care professionals. Patient circumstances will vary and some information may have become outdated as a result of more recent medical developments.

Ingredients	Mfr	Lot #	Expiry Date	Quantity	Measured	Checked
Clonazepam 2 mg Tablets	HLR			5		
Ora-Plus®	Paddock			50mL		
Ora-Sweet®*	Paddock			50mL		
OraPlus®/OraSweet® Mixture	SickKids			qs 100mL		

Can substitute Ora Sweet SF® for OraSweet®

### Equipment

- Mortar and pestle
- 125mL graduated measure x2
- Glass stirring rod
- 100mL clean amber bottle x1
- Shrink seal


### ✓ Procedure

1. Measure out and mix together Ora-Plus® and Ora-Sweet®
2. Crush tablets in the mortar with pestle to a fine powder. Levigate powder with small amounts Ora-Plus®/Ora-Sweet® vehicle until a smooth paste is formed. Add small amounts of vehicle to the paste until a liquid is formed. Transfer the contents into the graduate. Use additional vehicle to rinse the remaining drug from the mortar and add to the graduate.
3. Make up to the final volume of 100mL with the vehicle. Stir well.
4. Transfer the suspension into the amber glass bottle.
5. Label and affix shrink seal.

Storage: Room temperature

Expiry: 60 days

Sample Label:

	<b>CLONAZEPAM 0.1MG/ML ORAL SUSPENSION</b>
L#:	EXP:
ROOM TEMPERATURE	SHAKE WELL (100mL)

Date and Signature (Technician and Pharmacist Final Check): \_\_\_\_\_

\_\_\_\_\_