

# Dipyridamole

## 10mg/mL Oral Suspension

Last updated: April 2007

**Batch No:** \_\_\_\_\_

Disclaimers: CAUTION: The information on The Hospital for Sick Children World Wide Web Site has been developed specifically for The Hospital for Sick Children. It is provided for informational purposes for qualified health care professionals. This material is not intended as a substitute for consulting qualified health care professionals. Patient circumstances will vary and some information may have become outdated as a result of more recent medical developments.

Ingredients	Mfr	Lot #	Expiry Date	Quantity	Measured	Checked
Dipyridamole 50mg Tablets				20		
Ora-Plus®	Paddock			50mL		
*Ora-Sweet®	Paddock			50mL		
Ora-Plus®/Ora-Sweet®	SickKids			qs		
Mixture				100mL		

\*Can substitute with Ora-Sweet SF® for a sugar-free suspension.

### Equipment

- Glass mortar and pestle
- 125mL graduated measure x2
- Glass stirring rod
- 100mL clean amber bottle x1
- Shrink seal x1


### ✓ Procedure

1. Measure out and mix together the Ora-Plus® and Ora-Sweet® vehicles.
2. Put the Dipyridamole tablets into the mortar and let soak with a small amount of Ora-Plus®/Ora-Sweet® mixture until the coating on the tablets has softened.
3. Once softened, crush the tablets into a smooth paste. Add small portions of vehicle to the paste until a liquid is formed. Transfer the contents of the mortar to the graduate. Use additional vehicle to rinse the remaining drug from the mortar.
4. Make up to the final volume with vehicle. Stir well.
5. Transfer the suspension to the amber bottle.
6. Label and affix shrink seal.

Storage: Room temperature

Expiry: 60 days

Sample Label:

	<b>DIPYRIDAMOLE 10MG/ML ORAL SUSPENSION</b>
L#:	EXP:
ROOM TEMPERATURE	SHAKE WELL (100mL)

Date and Signature (Technician and Pharmacist Final Check): \_\_\_\_\_  
 \_\_\_\_\_