

Losartan

2.5mg/mL Oral Suspension

Last updated: April 2007

Batch No: _____

Disclaimers: CAUTION: The information on The Hospital for Sick Children World Wide Web Site has been developed specifically for The Hospital for Sick Children. It is provided for informational purposes for qualified health care professionals. This material is not intended as a substitute for consulting qualified health care professionals. Patient circumstances will vary and some information may have become outdated as a result of more recent medical developments.

| Ingredients | Mfr | Lot # | Expiry Date | Quantity | Measured | Checked |
|---------------------------------|---------------|-------|-------------|----------|----------|---------|
| Losartan 50mg Tablets (Cozaar®) | Merck Frosst* | | | 5 | | |
| Sterile Water for Irrigation | | | | 5mL | | |
| Ora-Plus® | Paddock | | | 50mL | | |
| Ora-Sweet® | Paddock | | | 50mL | | |
| OraPlus/OraSweet mixture | SickKids | | | qs 100mL | | |

*no substitution

Equipment

- Mortar and pestle
- 5mL oral syringe
- Stirring rod
- 2x 100mL graduate
- Amber bottle


✓ Procedure

1. Measure out and mix together the Ora-Plus® and Ora-Sweet®.
2. Count out Losartan 50mg tablets and place in a mortar.
3. Add water and let tablets soak.
4. Crush and levigate tablets to make a smooth paste.
5. Add vehicle in small amounts to the paste, mixing well until a liquid is formed. Transfer the contents to the graduate. Use additional vehicle to rinse the remaining drug from the mortar and add to the graduate
6. Make up to the final volume in the graduate using more vehicle. Stir well. (You will have a bit of the vehicle left over).
7. Transfer suspension to an amber bottle.
8. Label

Storage: Refrigerate

Expiry: 28 days

Sample Label:

| | |
|---|--|
|  | LOSARTAN 2.5MG/ML ORAL SUSPENSION |
| L# : | EXP: |
| REFRIGERATE | SHAKE WELL (100ML) |

Date and Signature (Technician and Pharmacist Final Check): _____
