



# Metronidazole 15mg/mL Oral Suspension

Last updated: November 2008

Batch No: \_\_\_\_\_

Disclaimers: CAUTION: The information on The Hospital for Sick Children World Wide Web Site has been developed specifically for The Hospital for Sick Children. It is provided for informational purposes for qualified health care professionals. This material is not intended as a substitute for consulting qualified health care professionals. Patient circumstances will vary and some information may have become outdated as a result of more recent medical developments.

Ingredients	Mfr	Lot #	Expiry Date	Quantity	Measured	Checked
Metronidazole 250 mg tablets	APOTEX			15		
Sterile Water for Irrigation USP				5 mL		
Chocolate Syrup	SickKids			qs 250mL		

## Equipment

- small mortar and pestle
- glass stirring rod
- 250 mL graduated measure
- 250 mL clean amber bottle x4
- shrink seal

## ✓ Procedure

1. Place tablets and water in the mortar and let sit for a few seconds. Levigate into a smooth paste using Chocolate Syrup.
2. Add more Chocolate Syrup to the paste until a liquid is formed. Transfer the contents to the graduated measure. Use additional Chocolate Syrup to rinse remaining drug from the mortar.
3. Make up to final volume with Chocolate Syrup. Stir well.
4. Transfer the suspension to the amber bottle.
5. Label and affix shrink seal.

Storage: Refrigerate

Expiry: 60 Days

	<b>METRONIDAZOLE 15MG/ML ORAL SUSPENSION</b>
L#:	EXP:
REFRIGERATE	SHAKE WELL (250mL)

Date and Signature (Technician and Pharmacist Final Check): \_\_\_\_\_  
\_\_\_\_\_