



# Midazolam HCl 3mg/mL Oral Syrup

Last updated: November 2008

Batch No: \_\_\_\_\_

Disclaimers: CAUTION: The information on The Hospital for Sick Children World Wide Web Site has been developed specifically for The Hospital for Sick Children. It is provided for informational purposes for qualified health care professionals. This material is not intended as a substitute for consulting qualified health care professionals. Patient circumstances will vary and some information may have become outdated as a result of more recent medical developments.

Ingredients	Mfr	Lot #	Expiry Date	Quantity	Measured	Checked
Midazolam HCl 5 mg/mL injection	HLR or Sandoz			300mL		
Chocolate Syrup	SickKids			qs 500mL		

## Equipment


- 500mL graduated measure
- 60mL injectable syringe x2
- 18g needle x2
- 25mL clean amber bottle x10
- 250mL clean amber bottle
- shrink seal x11

## ✓ Procedure

- 1. Withdraw 300mL of Midazolam, using the 60mL syringes and the 18g needles. Transfer the injection into the graduated measure.
- 2. Make up to the final volume of 500mL with Chocolate Syrup. Stir well.
- 3. Transfer the contents of the graduate into the amber bottles.
- 4. Label and affix shrink seal.
- 5. Label the remainder of the midazolam injection with a 30 day expiry and "FOR ORAL USE ONLY".

Storage: Room Temperature

Expiry: 60 days

		<b>MIDAZOLAM HCL 3MG/ML ORAL SYRUP</b>	
L#:		EXP:	
ROOM TEMPERATURE		SHAKE WELL (500ML)	

Date and Signature (Technician and Pharmacist Final Check): \_\_\_\_\_  
\_\_\_\_\_